




SURGICAL MASKS

Surgical masks are mainly used in the medical sector, such as doctors' surgeries, clinics or in nursing. Wearing them in everyday life can also fulfil official regulations. The mask is made of non-woven fabric and filter cloth. It has a flexible nose clip and an ear hook. Wearing the masks area-wide can prevent the spread of saliva or breath droplets of the wearer and thus the spread of the coronavirus.



01

Nonwoven
prevents droplets
from entering the
mask

02

Meltblown
barrier rate of
5 µm, particles
> 90 %

03

Nonwoven
absorbs moisture
close to mouth and
nose

Skin care fabric
Soft against skin,
breathable and light



SURGICAL MASKS

non-sterile

PERFORMANCE	STANDARD	VALUE
Nose clip length	YY 04069-2011	≥ 8 cm
Breaking strength of mask belt connection	YY 04069-2011	≥ 10 N
BFE	EN 14683-2019	≥ 98 %
PFE	YY 04069-2011	≥ 30 %
Respiratory resistance	EN 14683-2019	< 60 Pa
Synthetic blood penetration resistance	ISO 22609:2009	120 mmHg
Microbiological index	ISO 11737-1:2008	≤ 30 cfu/g
Biocompatibility	EN ISO 10993-1:2009	Non-toxic and non-irritating to human skin
Flame retardant performance	YY 04069-2011	Burning time less than 5 sec
Mask performance level	EN 0469-2011	Type IIR

PRODUCT INFORMATION

Dimensions	17,5 x 9,5 cm
Layers	3
Shelf life	36 months
Package	50 pcs per box
Storage	It is recommended to store in a dry, ventilated and non-corrosive gas environment



DISTRIBUTED BY:

hawo GmbH
Obere Au 2-4 | 74847 Obrigheim / Germany
T +49 (0) 6261 / 9770 - 0
F +49 (0) 6261 / 9770 - 69
info@hawo.com
www.hawo.com

Made in China

Technical modifications reserved
04/2020 | 9.344.101
Version 1.0